BOSCH

vivalytic

SARS-CoV-2 Test
*Severe Acute Respiratory Syndrome
Coronavirus 2*

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Package Contents

15 Vivalytic SARS-CoV-2 test cartridges for the detection of SARS-CoV-2 specific E gene sequence.

Safety Information

These instructions for use contain test-specific information only. For additional warnings and instructions, refer to the instructions for use provided with your Vivalytic analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the analyser. Make sure to insert samples properly. Remove spilled samples with a cloth soaked with the appropriate solution. You can use 70 % ethanol. If you need to disinfect the cartridge by spraying (for example after heavy spillage of sample), do not use the cartridge afterwards since it might be damaged by disinfectant.

WARNING

- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not use sample types, media and volumes that are not applicable for the test.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Handle potentially infectious samples according to laboratory standards and dispose samples according to regional and laboratory standards. Make sure to wear appropriate personal protective equipment (PPE).
- Place the cartridge on a clean surface.
- Be compliant with the national safety regulations and practices.

Additional Equipment (required but not provided)

Bosch Vivalytic one analyser Pipette (1000 µl)

Intended Use

The Bosch Vivalytic SARS-CoV-2 test is a qualitative PCR-based assay for the detection of SARS-CoV-2 from human nasopharyngeal and oropharyngeal swabs to aid in the diagnosis of a SARS-CoV-2 infection. Intended for use with a Vivalytic analyser by healthcare professionals.

Test Principle

Vivalytic SARS-CoV-2 is a real-time based PCR test utilizing reverse transcription (RT).

Storage und Usage Conditions

Storage and usage conditions can be taken either from the cartridge, pouch or box label.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, reverse transcription, DNA amplification and detection.

Reagents are bead, binding buffer, washing buffer and elution buffer.

Sample Type/Medium

The test is intended for use with nasopharyngeal and oropharyngeal swab samples in eNAT™ medium (COPAN Italia S.P.A.). Collect and store samples as indicated in the manufacturer's data sheet.

Sample Preparation

Shake the sample tube containing the swab sample and the liquid medium for homogenization and fill 300 µl of homogenized patient sample in the sample input of the cartridge. Do not use viscous samples that are difficult to pipette.

Vivalytic SARS-CoV-2 Test – Instructions for use

Test Result

The test result is shown on the screen.

Detection of the human control shows a successful extraction procedure and excludes an inhibition of the PCR reaction.

A test is rated as invalid only if both criteria are fulfilled; no detection of target RNA and no detection of human control.

In the printed test report, pathogen, results, control and information on user, patient and analyser are listed with a signature field. In case of positive results, the respective C_q value is displayed on the report.

In case of an invalid test, check if any notices are displayed after the run. Possible reasons for an invalid run might be poor sample quality or no RNA. Repeat the analysis with a new aliquot of the same sample if required.

Pay attention to correct sample type, the right sample collection and storage of the sample and cartridges prior to the test run.

In case of a failed test, first check for correct operating conditions of the analyser (refer to analyser's instructions for use, chapters' device safety information and technical data). Restart the analyser. If the problem persists, contact the customer service.

Quality Control

If required by your local or laboratory standards, quality control testing has to be performed. You can either use pre-characterized patient samples that were investigated by a reference testing method or purchase quality control materials. In case of unexpected results, repeat the analysis with another sample. If the result of a quality negative control sample remains positive, the analyser or its environment might be contaminated. Stop using the analyser and call the customer service. In case of repeated negative results for positive quality control samples, call also the customer service.

Limitations

The results of the Vivalytic SARS-CoV-2 test should be interpreted by a trained healthcare professional. The results of the Vivalytic SARS-CoV-2 test should not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogen being present in the sample at a level below assay sensitivity or a pathogen that is not covered by this assay.
- There is a risk of false negative values resulting from improperly collected, transported, or handled samples.

Analytical Sensitivity (Limit of Detection, 95 % Detection rate)

The concentration at a detection rate of 95 % (LoD) was determined (table 1).

Inclusivity and Exclusivity

Specificity was ensured by the selection of primers and their *in silico* analysis for possible cross-reactions based on publicly available nucleic acid sequences derived from the NCBI database.

To evaluate inclusivity, positive SARS-CoV-2 material (table 2) were purified and spiked into a PCR reaction and processed by a workflow using liquid components.

To exclude cross-reactivity (exclusivity), various strains of microorganisms representing common respiratory pathogens or closely related species were tested (table 3) by a workflow using liquid components. There was no evidence of microbial interference.

Interferences

Interferences were evaluated for endogenous and exogenous substances (table 4) that are potentially present in the patient sample. No interferences were detected.

Clinical Sensitivity and Specificity

The results derived from patient samples (positive and negative samples in Copan eNAT™ medium) collected in a clinical setting were compared with those of a reference method (table 5). In total, 50 samples were analysed.

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Bosch Healthcare Solutions GmbH